SKINNERS[°] 2.0 Compression

K

Sizing Guide

- 1. Print to Actual size / 100% scale Paper size: A4 (210×297mm)
- 2. Check the printed scale against the scale confirmation line below with a ruler.
- 3. Lay the paper on a flat surface and stand with your heel on the bottom line.
- 4. Mark the line above the longest toe.
- 5. Choose your size.
- ! Between sizes? Customers prefer to size up. In doubts? Leave us a message via e-mail or social media.

ALIGN YOUR HEEL

5cm -

SCALE CONFIRMATION When printed to 100% scale

47.₄₉

36-3;

Wiath Due to the high Rexibility and adaptability in rormation is not adaptability and adap

ty of the product, this information is not as suitable for customers with narrower feet. ^{essential} as for shoes. Suitable for customers with narrower feet. E.g.: Skinners size L is suitable for II cm wide 9.: Skinners Size L IS SUITAble for Lich Vider toe box).

3ICM 12 1/5,

^{29.5}cm ^[] 2/3,

28cm

17.

27₀m 10 5/8"

^{25.5}cm In.

4cm 1/2

the line above must scale 5cm

the line above must scale 2" When printed to 100% scale SCALE CONFIRMATION

5"

EU	US		Foot Length	
Size	Women	Men	cm	inch
36-37	5.5-6.5	4.5-5.5	22.5-24	81/8-91/2
38-39	7–8	6–7	24-25.5	9½-10
40-42	8.5-10.5	7.5-9.5	25.5–27	10-105%
43-44	11-12	10-11	27–28	105%-11
45-46	12.5-13.5	11.5-12.5	28-29.5	11-11²⁄3
47-49	14-15.5	13-14.5	29.5-31	11 ² / ₃ -12 ¹ / ₅

E-mail: hello@skinners.cc FB: skinnersoriginal IG: skinners_footwear